


Ambler is excited to introduce our newest product the “Perfect Pupil”

The Flexible Iris Expansion Ring provides your surgeon with the best pupil dilation to treat intraoperative floppy iris syndrome (IFIS) patients more proficiently.

- Easy to insert using the Perfect Pupil Injector
- Securely expands and holds a billowing floppy iris through 315°
- Keeps pupil dilated to 7mm even during progressive miosis
- Provides full control during rhexis creation
- Aids in cortex and epithelial cell removal
- Retains pupil’s normal size, shape and function postoperatively



To Order Call Toll Free: 1-888-407-0006 Today!

The Perfect Pupil Injectable is a flexible iris expansion ring with a 0.24mm groove and 9 scalloped tabs to securely hold a billowing, floppy iris through 315 degrees, made of flexible polyurethane to keep a small and constricting pupil dilated to 7mm throughout the surgical procedure, CE marked, FDA approved, sold per each, sterile, disposable.



Compare the “Perfect Pupil” to the Competition:

	Perfect Pupil	Malyugin Ring*	Morcher Pupil Dilator*
Provided Sterile	Yes	No	Yes
Dilation Size	7.0mm	6.00-6.50mm	5.00-6.00mm
Easy Unloading Arm	Yes	No	No
Smallest Incision Size	2.75mm	2.20mm	2.85mm
Pupil Ring Injector	Optional	Disposable	Optional
Injector Price	\$495.00	Included	\$705.00
Price Each	\$95.00	\$125.00	\$149.00
Price 6 / Pack	\$475.00	\$750.00	\$894.00

*The above referenced companies/products are registered trademarks of the respective owners. Ambler Surgical has no affiliation with these companies, nor do we distribute their products.

Ambler Surgical Toll Free: 1-888-407-0006
Order Online www.AmblerSurgical.com